

ALLERGUARD^{ES}

NU-608 CLASS I ANIMAL BEDDING DISPOSAL CABINET



- Personnel & Environmental Protection
- Maintains an Active Air Barrier with a Constant 105 FPM [0.53 m/s] Inflow Air Velocity
- Easy-to-Clean Stainless Steel Work Surface
- Pneumatic Chute Lift for Easy Bag Change Access
- Stainless Steel Bang-Bar Over the Trash Chute
- Deodorizing Charcoal Exhaust Filter
- Heavy-Duty Reinforced Casters



Easily attach garbage bag to waste chute. No tools necessary!



Extra large 14 in [356 mm] access opening comfortably fits rat cages.




Stainless steel bang bar ensures waste gets successfully dumped into chute.

ALLEREGARD^{ES}

NU-608 CLASS I ANIMAL BEDDING DISPOSAL CABINET

AeroMax™ Electronic Control System



NU-608 FEATURES & SPECIFICATIONS

Class I Animal Bedding Disposal Cabinet

Model #	NU-608-400
Cabinet Style	Animal Bedding Disposal Station (Console)
Protection	Personnel, Environment
Nominal Width	4 ft [1.2 m]
Workspace Dimensions (W x D x H)	46.25 x 29.75 x 28 in [1174 x 756 x 711 mm]
Exterior Dimensions (W x D x H)	52.75 x 33.375 x (84.75 - 96.75) in [1340 x 848 x (2153 - 2457) mm]
Net Weight	400 lbs [181 kg]
Certifications	UL listed, ADA Compliant
Access Opening	14 in [356 mm]
Electrical Requirements	115 VAC / 60 Hz
Air Barrier	105 fpm [0.53 m/s] Inflow
Filter Type	6 in [152 mm] Single HEPA 99.99% @ 0.3 micron Efficiency
Exhaust Module	C-25 Charcoal Exhaust Filter Module
Control System	AeroMax™ w/ PressureFlow™ Monitor
Lighting	LED (50-150 fc [538 - 1600 lux])
Noise Level	Up to 65 dba
Work Zone Construction	Single Rectangular Disposal Chute w/ Stainless Steel Bang-Bar
Base Stand Options	Telescoping Base Stand w/ Casters
	Motorized Base Stand w/ Casters
	Standard Feature
	Optional Feature



To Learn More or Request a Quote:

Visit: www2.nuaire.com/608-flyer

or Call: 1-763-553-1270

Toll Free: 1-800-328-3352

