







REV	ECO	DESCRIPTION	DATE	DFTM	CHKD
A	13949	RELEASED TO PRODUCTION	1/23/2019	TH	BK



MODEL DIMS			
MODEL	DIM "A"	DIM "B"	HEPA FILTER SIZE
NU-340-530	32 [813]	80 1/2 [2045]	58 [1473] X 30 [762] X 3 [76]
NU-340-536	38 [965]	86 1/2 [2197]	58 [1473] X 36 [914] X 3 [76]

\* OVERALL DIMENSIONAL TOLERANCE ±1/4" [6.35]  
ALL OTHER DIMENSIONS ±1/8" [3.17]

TITLE NU-340-500 HORIZONTAL AIRFLOW CONSOLE W/ CASTORS					
DFTM TH	DATE 1/23/2019	CHKD BK	SHEET 1	OF 1	
DRAWING NUMBER		BCD-20629		A	

SECTION A-A

INCH  
MILLIMETER