

V-1 plus, Personal Vortex

Basic Plus
Product Class

DESCRIPTION

Vortex **V-1 plus** is an ideal instrument for gentle mixing to vigorous resuspension of cells and biological and chemical liquid components in tubes using eccentric mechanism.

Vortex has two modes:

1. Continuous operation;
2. Impulse operation (activated by pressing the cap with the tube's bottom).

SPECIFICATIONS

Eccentric mixing principle	+
Speed control range	500-3000 RPM
Acceleration time	2 s
Maximum continuous operation time	24 h
Mixing module for tubes	from 0.2 to 50 ml
Maximum load	30 g
Orbit	4 mm
Overall dimensions (W×D×H)	90x150x80 mm
Weight	0.8 kg
Input current/power consumption	12 V, 320 mA / 3.8 W
External power supply	Input AC 100–240 V; 50/60 Hz; Output DC 12 V



CAT. NUMBER

BS-010203-AAG	230VAC 50/60Hz Euro plug
BS-010203-AAK	100-240VAC 50/60Hz Multi plug (EU, UK, AU, US)
BS-010203-BK	IQ OQ document
BS-010203-CK	PQ document